

PIPEJET® NANODISPENSER KIT



PIPEJET NANODISPENSER KIT

The PipeJet Nanodispenser Kit includes everything you need to realize your application. It includes a PipeJet dispensing module and the control electronics as well as the intuitive PipeJet control software.

In order to allow the simplest possible integration of the PipeJet Nanodispenser, the control electronics processes external trigger signals. Additionally, the complete PipeJet functionality can be embedded into your system via the provided DLLs.

Applications

Aliquoting & distributing precious samples (Enzymes, Proteins etc.)

Dispensing of cells
Dispensing of aggressive media which degrades tubing and valves of conventional system

Aspiration from and pipetting into 1536 well-MTP

Materials

All liquid contaminated materials:

- are low cost disposables
- can be sterilized
- are USP class VI compliant

Customized materials and dispenser modules are possible



BioFluidix GmbH
Engesserstr. 4a · 79108 Freiburg
Fon: +49 761 458938 - 0 · Fax: +49 761 458938 - 99
info@biofluidix.com · www.biofluidix.com

